

- e-CHECK*** (XS) Quality Control for the **ProCyte Dx*** Analyzer
- e-CHECK*** (XS) contrôle qualité pour l'analyseur **ProCyte Dx***
- e-CHECK*** (XS) Control de calidad para el analizador **ProCyte Dx***
- e-CHECK*** (XS) Qualitätskontrolle für das **ProCyte Dx*** Analysegerät
- e-CHECK*** (XS) Controllo qualità per l'analizzatore **ProCyte Dx***
- e-CHECK*** (XS) IDEXXプロサイトDx*の精度管理

L1 (Low / Bas / Bajo / Niedrig / Basso / 低濃度)

LOT



Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)		
HCT (%)		
HGB (g/dL)		
MCV/VGM (fL)		
MCH/TCMH (pg)		
MCHC/CCMH (g/dL)		
PLT (K/μL)		
RDW-SD (fL)		
RDW-CV (%)		
PDW (fL)		
PCT (%)		
MPV (fL)		
P-LCR (%)		
WBC/GB (K/μL)		
NEUT# (K/μL)		
LYMPH# (K/μL)		
MONO# (K ⁹ /μL)		
EO# (K/μL)		
RBC-O (10 ⁶ /μL)		
PLT-O (K/μL)		
RETIC (K/μL)		
% RETIC (%)		

L2 (Normal / Normale / 正常濃度)

LOT

52730805



12/21/25

Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)	4.36	4.09 - 4.63
HCT (%)	36.4	32.9 - 39.9
HGB (g/dL)	12.1	11.6 - 12.6
MCV/VGM (fL)	83.6	75.6 - 91.6
MCH/TCMH (pg)	27.8	24.1 - 31.5
MCHC/CCMH (g/dL)	33.2	27.1 - 39.3
PLT (K/μL)	205	160 - 250
RDW-SD (fL)	44.9	37.9 - 51.9
RDW-CV (%)	15.9	12.9 - 18.9
PDW (fL)	9.5	6.5 - 12.5
PCT (%)	0.19	0.11 - 0.27
MPV (fL)	9.4	6.4 - 12.4
P-LCR (%)	14.0	8.1 - 19.9
WBC/GB (K/μL)	7.23	6.48 - 7.98
NEUT# (K/μL)	3.49	2.43 - 4.55
LYMPH# (K/μL)	2.30	1.68 - 2.92
MONO# (K ⁹ /μL)	0.69	0.00 - 1.45
EO# (K/μL)	0.75	0.39 - 1.11
RBC-O (10 ⁶ /μL)	4.17	3.47 - 4.87
PLT-O (K/μL)	203	153 - 253
RETIC (K/μL)	120.0	90.0 - 150.0
% RETIC (%)	2.75	1.95 - 3.55