

- e-CHECK*** (XS) Quality Control for the **ProCyte Dx*** Analyzer
- e-CHECK*** (XS) contrôle qualité pour l'analyseur **ProCyte Dx***
- e-CHECK*** (XS) Control de calidad para el analizador **ProCyte Dx***
- e-CHECK*** (XS) Qualitätskontrolle für das **ProCyte Dx*** Analysegerät
- e-CHECK*** (XS) Controllo qualità per l'analizzatore **ProCyte Dx***
- e-CHECK*** (XS) IDEXXプロサイトDx*の精度管理

L1 (Low / Bas / Bajo / Niedrig / Basso / 低濃度)

LOT



Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)		
HCT (%)		
HGB (g/dL)		
MCV/VGM (fL)		
MCH/TCMH (pg)		
MCHC/CCMH (g/dL)		
PLT (K/μL)		
RDW-SD (fL)		
RDW-CV (%)		
PDW (fL)		
PCT (%)		
MPV (fL)		
P-LCR (%)		
WBC/GB (K/μL)		
NEUT# (K/μL)		
LYMPH# (K/μL)		
MONO# (K ⁹ /μL)		
EO# (K/μL)		
RBC-O (10 ⁶ /μL)		
PLT-O (K/μL)		
RETIC (K/μL)		
% RETIC (%)		

L2 (Normal / Normale / 正常濃度)

LOT

60200805



04/12/26

Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)	4.34	4.07 - 4.61
HCT (%)	35.2	31.7 - 38.7
HGB (g/dL)	11.7	11.2 - 12.2
MCV/VGM (fL)	81.0	73.0 - 89.0
MCH/TCMH (pg)	27.0	23.3 - 30.7
MCHC/CCMH (g/dL)	33.2	27.1 - 39.3
PLT (K/μL)	208	163 - 253
RDW-SD (fL)	55.9	48.9 - 62.9
RDW-CV (%)	20.4	17.4 - 23.4
PDW (fL)	10.1	7.1 - 13.1
PCT (%)	0.19	0.11 - 0.27
MPV (fL)	9.4	6.4 - 12.4
P-LCR (%)	15.3	9.4 - 21.2
WBC/GB (K/μL)	7.30	6.55 - 8.05
NEUT# (K/μL)	3.63	2.57 - 4.69
LYMPH# (K/μL)	2.10	1.48 - 2.72
MONO# (K ⁹ /μL)	0.79	0.03 - 1.55
EO# (K/μL)	0.78	0.42 - 1.14
RBC-O (10 ⁶ /μL)	4.23	3.53 - 4.93
PLT-O (K/μL)	166	116 - 216
RETIC (K/μL)	126.4	96.4 - 156.4
% RETIC (%)	2.91	2.11 - 3.71