

- e-CHECK*** (XS) Quality Control for the **ProCyte Dx*** Analyzer
- e-CHECK*** (XS) contrôle qualité pour l'analyseur **ProCyte Dx***
- e-CHECK*** (XS) Control de calidad para el analizador **ProCyte Dx***
- e-CHECK*** (XS) Qualitätskontrolle für das **ProCyte Dx*** Analysegerät
- e-CHECK*** (XS) Controllo qualità per l'analizzatore **ProCyte Dx***
- e-CHECK*** (XS) IDEXXプロサイトDx*の精度管理

L1 (Low / Bas / Bajo / Niedrig / Basso / 低濃度)

LOT



Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)		
HCT (%)		
HGB (g/dL)		
MCV/VGM (fL)		
MCH/TCMH (pg)		
MCHC/CCMH (g/dL)		
PLT (K/μL)		
RDW-SD (fL)		
RDW-CV (%)		
PDW (fL)		
PCT (%)		
MPV (fL)		
P-LCR (%)		
WBC/GB (K/μL)		
NEUT# (K/μL)		
LYMPH# (K/μL)		
MONO# (K ⁹ /μL)		
EO# (K/μL)		
RBC-O (10 ⁶ /μL)		
PLT-O (K/μL)		
RETIC (K/μL)		
% RETIC (%)		

L2 (Normal / Normale / 正常濃度)

LOT

62150805



10/25/26

Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)	4.31	4.04 - 4.58
HCT (%)	35.1	31.6 - 38.6
HGB (g/dL)	11.7	11.2 - 12.2
MCV/VGM (fL)	81.5	73.5 - 89.5
MCH/TCMH (pg)	27.1	23.4 - 30.8
MCHC/CCMH (g/dL)	33.3	27.2 - 39.4
PLT (K/μL)	210	165 - 255
RDW-SD (fL)	46.6	39.6 - 53.6
RDW-CV (%)	17.0	14.0 - 20.0
PDW (fL)	9.7	6.7 - 12.7
PCT (%)	0.19	0.11 - 0.27
MPV (fL)	9.2	6.2 - 12.2
P-LCR (%)	13.8	7.9 - 19.7
WBC/GB (K/μL)	7.07	6.32 - 7.82
NEUT# (K/μL)	3.47	2.41 - 4.53
LYMPH# (K/μL)	2.15	1.53 - 2.77
MONO# (K ⁹ /μL)	0.70	0.00 - 1.46
EO# (K/μL)	0.75	0.39 - 1.11
RBC-O (10 ⁶ /μL)	4.24	3.54 - 4.94
PLT-O (K/μL)	171	121 - 221
RETIC (K/μL)	114.2	84.2 - 144.2
% RETIC (%)	2.65	1.85 - 3.45