

- e-CHECK*** (XS) Quality Control for the **ProCyte Dx*** Analyzer
- e-CHECK*** (XS) contrôle qualité pour l'analyseur **ProCyte Dx***
- e-CHECK*** (XS) Control de calidad para el analizador **ProCyte Dx***
- e-CHECK*** (XS) Qualitätskontrolle für das **ProCyte Dx*** Analysegerät
- e-CHECK*** (XS) Controllo qualità per l'analizzatore **ProCyte Dx***
- e-CHECK*** (XS) IDEXXプロサイトDx*の精度管理

L1 (Low / Bas / Bajo / Niedrig / Basso / 低濃度)

LOT



Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)		
HCT (%)		
HGB (g/dL)		
MCV/VGM (fL)		
MCH/TCMH (pg)		
MCHC/CCMH (g/dL)		
PLT (K/μL)		
RDW-SD (fL)		
RDW-CV (%)		
PDW (fL)		
PCT (%)		
MPV (fL)		
P-LCR (%)		
WBC/GB (K/μL)		
NEUT# (K/μL)		
LYMPH# (K/μL)		
MONO# (K ⁹ /μL)		
EO# (K/μL)		
RBC-O (10 ⁶ /μL)		
PLT-O (K/μL)		
RETIC (K/μL)		
% RETIC (%)		

L2 (Normal / Normale / 正常濃度)

LOT

62710805



12/20/26

Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)	4.33	4.06 - 4.60
HCT (%)	35.8	32.3 - 39.3
HGB (g/dL)	12.1	11.6 - 12.6
MCV/VGM (fL)	82.7	74.7 - 90.7
MCH/TCMH (pg)	27.9	24.2 - 31.6
MCHC/CCMH (g/dL)	33.8	27.7 - 39.9
PLT (K/μL)	209	164 - 254
RDW-SD (fL)	44.3	37.3 - 51.3
RDW-CV (%)	15.9	12.9 - 18.9
PDW (fL)	8.5	5.5 - 11.5
PCT (%)	0.17	0.09 - 0.25
MPV (fL)	8.3	5.3 - 11.3
P-LCR (%)	8.3	2.4 - 14.2
WBC/GB (K/μL)	7.01	6.26 - 7.76
NEUT# (K/μL)	3.37	2.31 - 4.43
LYMPH# (K/μL)	2.25	1.63 - 2.87
MONO# (K ⁹ /μL)	0.66	0.00 - 1.42
EO# (K/μL)	0.73	0.37 - 1.09
RBC-O (10 ⁶ /μL)	4.29	3.59 - 4.99
PLT-O (K/μL)	175	125 - 225
RETIC (K/μL)	112.5	82.5 - 142.5
% RETIC (%)	2.60	1.80 - 3.40